Serial No.: 10/568,544

IN THE CLAIMS:

Please amend the claims as follows:

- 1. (Currently amended) A method for examining colorectal cancer and colorectal adenoma comprising assaying the reactivity <u>binding</u> of an antibody to <u>an antigen in</u> tissues, body fluid or feces of patients, or extracts thereof, wherein said antibody-reacts <u>with an antigen</u> [[that]] is present in cells expressing HEC-GlcNAc6ST gene encoding GlcNAc-6-sulfotransferase and [[that]] is absent <u>or almost absent</u> in cells expressing GlcNAc6ST-1 or I-GlcNAc6ST gene.
- (Currently amended) The method of claim 1, wherein said antigen is present in cells transduced with HEC-GlcNAc6ST gene and is absent er-almost absent-in cells transduced with GlcNAc6ST-1 gene or I-GlcNAc6ST gene.
- 3. (Currently amended) The method of claim 1 [[or 2]], wherein said antigen comprises the sugar residues expressed by the following formula:

R1-Gal β1-3/4(SO₃-6) GlcNAc β1-R2

where, R1 represents sugar residues added by other enzymes and is not limited in structure[[.]]; Gal β represents β galactose, GlcNac β represents β N-acetylglucosamine, Gal β 1-3/4 represents binding of 1 position of Gal β and 3 position and/or 4 position of GlcNac β , (SO₃-6) represents addition of a sulfate group to 6 position of GlcNac β , R2 represents –3GalNac α , -3Gal β or –2Man α and binds to 1 position of GlcNac β .

- 4. (Currently amended) The method <u>according to claim 1 of any one of claims 1 to 3</u>, wherein said antibody is MECA-79 antibody (Pharmingen, catalog No. 09961D).
- (Currently amended) A method for examining colorectal cancer and colorectal adenoma comprising assaying the reactivity <u>binding</u> of MECA-79 antibody or <u>another</u>

Serial No.: 10/568,544

antibody having the same binding specificity its equivalent with tissues, body fluid, feces or extract thereof of test subjects.

- 6. (Currently amended) The method <u>according to claim 5</u> of any one of claims 1 to 5 comprising reacting <u>binding</u> a labeled probe to said antibody and assaying the label qualitatively or quantitatively.
- 7. (Withdrawn) An examination reagent for colorectal cancer and colorectal adenoma comprising, as a major component, an antibody (including MECA-79 antibody) reacting specifically with an antigen carrying sugar residues, which is present in cells expressing HEC-GlcNAc6ST gene and is absent or almost absent in cells expressing GlcNAc6ST-1 or GlcNAc6ST gene.
- 8. (Withdrawn) An examination reagent for colorectal cancer and colorectal adenoma comprising, as a major component, an antibody (including MECA-79 antibody) reacting specifically with an antigen carrying sugar residues, which is present in cells transduced with HEC-GlcNAc6ST gene and is absent or almost absent in cells transduced with GlcNAc6ST-1 or GlcNAc6ST gene.
- 9. (Withdrawn) An examination reagent for colorectal cancer and colorectal adenoma comprising, as a major component, an antibody (including MECA-79 antibody) reacting specifically with an antigen carrying sugar residues, which are present in tissues, body fluid or feces of patients with colorectal cancer and colorectal adenoma and expressed by the following general formula:

R1-Gal β1-3/4 (SO₃-6) GlcNAc β1-R2

where, R1 represents sugar residues added by other enzymes and is not limited in structure[[.]; Gal β represents β galactose, GlcNAc β represents β N-acetylglucosamine, Gal β 1-3/4 represents binding of 1 position of Gal β and 3 position and/or 4 position of GlcNAc β , (SO₃-6) represents addition of a sulfate group to 6

Serial No.: 10/568,544

position of GlcNAc β , R2 represents –3GalNAc α , -3Gal β or –2Man α and binds to 1 position of GlcNAc β .

10. (Withdrawn) An examination reagent for colorectal cancer and colorectal adenoma comprising MECA-79 antibody as a major component.